



Standard Features:

- Compact, Hard-Mountable to Specialty Equipment
- 10,000 Footcandles
- Intense, Cold Illumination
- 4-Position, Solid State Intensity Control
- Light Weight & Rugged, High-Impact Plastic Housing
- 2 Year Warranty

Accessories & Options:

- Optional Single Gooseneck Light Guide
- Optional 60mm ID Ringlight
- Accepts SX-Series Light Guide Adapters

Model 3100 30 Watt Small-Footprint Illuminator


Technical Information

Input:	115 V AC, 60Hz
Lamp:	30 Watt, 10.8 V, Type EKZ qtz halogen
Color Temp.:	3100° Kelvin
Lamp Life:	200 hrs. full intensity, 10,000 at min. intensity
Lamp Output:	10,000 footcandles at fiber optic insertion plane
Fiber Optic Interface:	25mm (A Style)
Intensity Control:	4-Position Solid State Switch
Adapters	SX series adapter
Safety Certifications:	N/A
Dimensions:	3" x 4.5" x 6.5"; 4.5 lbs.

Options

3100 002831000000	30 Watt, 115 V AC Illuminator
3100-1 660000051001	30 Watt, 115 V AC Illuminator w/ BGG 1823M Single Gooseneck, (1) LH759 Focusing Lens
3100-A37P 660000051004	30 Watt, 115 V AC Illuminator w/ A3739P Ring Light, SX-10 Adapter

Model 3100 Dimensions:



Dolan-Jenner
industries

Dolan-Jenner Industries
159 Swanson Rd.
Boxborough, MA 01719 USA

Phone: 800.833.4237
Fax: 978.264.0292
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM