



Standard Features:

- Up To 2x The Light Output of Standard Quartz Halogen
- Brilliant White Light Provides True Color Rendering & Balance
- 100 Watt AC Regulated Power Supply
- 100 Watt Metal Halide Lamp
- 5500°K to 6500°K Color Temp.
- 0-100% Intensity Control Via Manual Iris
- External Access to Color Filter
- Integral UV/IR Filter on Visible Light Models
- Cool housing and nosepiece temp.
- UV Model w/o Filter for Use w/ Quartz Fiber Optics for Use in Some UV Applications
- UL/C-UL Listed
- 2 Year Warranty

Accessories & Options:

- 35mm Color Filters
- External Filter Holder Standard
- Manually Adjustable Iris Standard

Model MH100 Metal Halide Fiber Optic Illuminator

Technical Information



Input:	115V AC 50/60 Hz; 230V AC 50/60 Hz
Lamp:	100 Watt Metal Halide
Lamp Output:	Up to 3000 Lumens at fiber optic insertion plane
Color Temperature:	6000° Kelvin at max. intensity
Lamp Life:	3000 Hours (Nominal to 60% initial intensity)
Fiber Optic Interface:	MH100A - 0.985 in. dia.
Light Output Regulation:	AC voltage regulated
Intensity Control:	Manual Iris 0-100%
Dimensions:	7.25" W x 9.50" D x 4.60" H; 12 lbs

Model & P/N

Option Description

MH100A 660001050102	115V AC, 0.985" Nosepiece, Manual Iris, External Filter Holder, North American Line Cord
MH100A1 660001051102	230V AC, 0.985" Nosepiece, Manual Iris, External Filter Holder, European Line Cord

MH100 Dimensions:



Dolan-Jenner
industries

Dolan-Jenner Industries
159 Swanson Rd.
Boxborough, MA 01719 USA

Phone: 800.833.4237
Fax: 978.264.0292
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM