



Standard Features:

- 10,000 Footcandles
- Intense, Cold Illumination
- Convection-Cooled
- Sleek, Stackable Housing
- 3-Year Warranty
- UL, C-UL, CE Approved

Accessories & Options:

- Optional 230 V AC Input
- Optional Single Gooseneck Light Guide
- Optional Filter Holder & IR Filter

Model MI-30 30 Watt Convection Cooled Illuminator





Technical Information

| | |
|------------------------|---------------------------------------|
| Input: | 115 V AC (230 V AC Optional) |
| Lamp: | 30 Watt, 10.8 V, Type EKZ qtz halogen |
| Color Temp.: | 3100° Kelvin |
| Lamp Life: | 200 hrs. full intensity |
| Lamp Output: | 10,000 Foot Candles |
| Fiber Optic Interface: | 0.59" |
| Noise Level: | Silent, No Fan, No Vibration |
| Safety Certifications: | UL, C-UL, CE |
| Dimensions: | 6.0" x 8.5" x 5.3"; 7.25 lbs. |

Options

| | |
|------------|---|
| MI-30: | 30 Watt, 115 V AC Illuminator |
| MI-301: | 30 Watt, 230 V AC Illuminator |
| MI-30F: | 30 Watt, 115 V AC Illuminator w/ Filter Holder & IR Filter |
| MI-301F: | 30 Watt, 230 V AC Illuminator w/ Filter Holder & IR Filter |
| MI-30SG: | 30 Watt, 115 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense |
| MI-30SG1: | 30 Watt, 230 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense |
| MI-30SGF: | 30 Watt, 115 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense, Filter Holder & IR Filter |
| MI-30SG1F: | 30 Watt, 230 V AC Illuminator w/ BGG1823M Single Gooseneck, (1) LH759 Focusing Lense, Filter Holder & IR Filter |

MI-30 Dimensions:



[mm]
in

Dolan-Jenner
industries

Dolan-Jenner Industries
159 Swanson Rd.
Boxborough, MA 01719 USA

Phone: 800.833.4237
Fax: 978.264.0292
Email: support@dolan-jenner.com

1.800.83.FIBER

WWW.DOLAN-JENNER.COM