

DOLAN-JENNER ADAPTERS LIGHT GUIDE TO LIGHT SOURCE

Model	Part Number	ID	OD	FOR "A" or "B" Style Nose	Note:
ST-1	SK38450001	0.191"	N/A	A	Fits model 190 only.
SX-5	020863031001	0.191"	0.985"	A	
SX-5B	010041060000	0.191"	0.587"	B	
D-5	002930060001	0.191" x 2	0.985"	A	Dual output adapter
SX5.5	020863181000	0.217"	0.985"	A	
SX-14	020863031005	0.250"	0.985"	A	
SX-6	020863031002	0.316"	0.985"	A	
SX-6A	010040020000	0.316"	0.587"	B	
SX-6B	010041070000	0.316"	0.587"	B	
D-6	002930060002	0.316" x 2	0.985"	A	Dual output adapter
SX-6.5	020863031004	0.381"	0.985"	A	
SX-7	020863031003	0.437"	0.985"	A	
SX-7BR	010040140002	0.440"	0.587"	B	Obsolete
SX-9	020863050001	0.506"	0.985"	A	
SX-10	020863050002	0.565"	0.985"	A	
SX-10BR	010040140000	0.565"	0.587"	B	Obsolete
SX-10.5	010040170000	0.590"	0.985"	A	Replaced by SX-13
SX-13	010041080000	0.590"	0.985"	A	Adapts Mi-150 tips to "A"
SX-11	020863050003	0.629"	0.985"	A	
SX-11BR	010040140001	0.629"	0.587"	B	Obsolete
SX-12	010040060005	0.722"/0.629"	0.985"	A	Accepts Schott/Fostec guides
SX-8	020863050004	0.758"	0.985"	A	
SX-8BR	010040140003	0.758"	0.587"	B	Obsolete

"A" NOSE = 1.00" ID, "B" NOSE = 0.590" ID