



MANUFACTURER	MANUFACTURER'S MODEL NUMBER	RINGLIGHT	ADAPTER SALES ORDER CODE
Aus Jena	M-1, M-2, M-3		360-10 MAP
Aus Jena	Technival 2/Citoval 2	A3739P, A3739M	NOT REQUIRED
LEICA Cambridge Instruments	(B&L) Stereo or Stereozoom 1 thru 5, GZ-4	A3739P, A3739M	360-4 MAP
LEICA Cambridge Instruments	(B&L) Stereozoom 6 & SZ-4	A3739P, A3739M	360-4A MAP
LEICA Cambridge Instruments	(B&L) Stereozoom 6 & SZ-4 Auxiliary Lenses .3X, .5X, .75X, 1.5X, 2.0X & Lens Shield	A3739P, A3739M	360-5 MAP
LEICA Cambridge Instruments	(B&L) Stereozoom 7	A3739P, A3739M	360-3 MAP
LEICA Cambridge Instruments	(B&L) Stereozoom 7, Low Magnification	A3739P, A3739M	360-3A MAP
LEICA Cambridge Instruments	All Cycloptic (AO) Series, Stereostar/Zoom Models w/.5X Auxiliary Lens	A3739P, A3739M	360-5 MAP
LEICA Cambridge Instruments	(AO) Stereostar/zoom Models w/1.25X, 1.5X or Lens Shield	A3739P, A3739M	360-6 MAP
LEICA Cambridge Instruments	(AO) Stereostar/zoom Models 569, 570, 580	A3739P, A3739M	360-6B MAP
LEICA Cambridge Instruments	GZ-6 .5X, 2.0X	LP25A36M	NOT REQUIRED
LEICA Cambridge Instruments	S4E, S6D/E, S8APO	LP25LI36M	NOT REQUIRED
LEICA Cambridge Instruments	M10	A3739P, A3739M	360-12A MAP
LEICA Cambridge Instruments	Leica S4/S6 ring light	LP25LI36M	NOT REQUIRED
Edmund Scientific	37659	A3739P, A3739M	360-6 MAP
Edmund Scientific	38334, 38335	A3739P, A3739M	360-6 MAP
Edmund Scientific	85056, 37691	A3739P, A3739M	360-12A MAP
Meiji	EMZ Series	A3739P, A3739M	NOT REQUIRED
NIKON	SMZ-U	700009570036	NOT REQUIRED
NIKON	SMZ10A	700009570036	NOT REQUIRED
NIKON	SM, SMZ, SMZ-10	A3739P, A3739M	360-5 MAP
NIKON	SMZ 645, SMZ660, SMZ660V	A645NK41	LOCKING RING INCLUDED
NIKON	SMZ1500	700009570036	NOT REQUIRED
NIKON	SMZ800	LPA25A36M	NOT REQUIRED
NIKON	SMZ1000 w/.5X to 1.5X Auxiliary Lenses	A645NK41	LOCKING RING INCLUDED
NIKON	SMZ1000 w/2X Auxiliary Lens	700009570036	NOT REQUIRED
NIKON	SMZ-1/1B, SMZ-DOD, SMZ-2B/2T	A3739P, A3739M	360-11 MAP
OLYMPUS	SZH	700009570036	NOT REQUIRED
OLYMPUS	SZH	A3739P, A3739M	NOT REQUIRED
OLYMPUS	SZ-II, SZ-III, SZ-III-Tr, JM/JM-Tr, SZ-30, SZ-40, SZ-60 Models with .5X & .75X Auxiliary Lenses	A3739P, A3739M	360-7 MAP
OLYMPUS	SZ-II, SZ-III, SZ-III-Tr, JM/JM-Tr W/O AUX LENS	A3739P, A3739M	360-6A MAP
OLYMPUS	SZ-30, SZ-40, SZ-60 Models without auxiliary lenses	AOLM48	LOCKING RING INCLUDED
OLYMPUS	SZX9, SZX12	700009570036	NOT REQUIRED
OLYMPUS	SZ-51, SZ-61	AOLM48	LOCKING RING INCLUDED
OLYMPUS	SZ-7	700009570036	NOT REQUIRED
OLYMPUS	VM, VMZ	A3739P, A3739M	360-12 MAP
POLAROID	Model MP4 Camera	A3739P, A3739M	360-7A MAP
UNITRON	FSB, ZSB & ZST	A3739P, A3739M	360-9 MAP
WILD	M-3, M3Z Plan, M-4, M-4A, M-5A, M-7A, M7S, M8,	A3739P, A3739M	NOT REQUIRED
WILD	M8, M400, MZB, Lower Magnification Only	A3739P, A3739M	360-12A MAP
WOLFE	Advanced & Zoom Models	A3739P, A3739M	360-9 MAP
ZEISS	SR Scopes & Aux Lenses: F=50, 100, 125, 175,	A3739P, A3739M	360-9 MAP
ZEISS	SV8 & IV/B	A3739P, A3739M	360-5 MAP
ZEISS	Microscopes with Objective O.D.s from 65-66.5mm	A3739P, A3739M	360-12A MAP
ZEISS	GSZ	A3739P, A3739M	NOT REQUIRED
ZEISS	STEMI 2000	A3739P, A3739M	NOT REQUIRED
UNITRON	ZST, ZSB 48" Ringlight	AUNI48	LOCKING RING INCLUDED

**150 WATT HALOGEN RING LIGHT ILLUMINATION SYSTEMS (115 & 230 volt input available)**

Scope Manuf.	MODEL	DESCRIPTION	WHERE USED	MOUNTING
NIKON	MI-150RA	MI-150 115V/RL 11 DJI SYS	fits Nikon SMZ645, 660	lock ring
UNIVERSAL	MI-150RB	MI-150 25 /RL SYS	66 mm ID universal fit, 25 deg transition	2 thumbscrews
UNIVERSAL	MI-150RC	MI-150 115V/RC 0 DJI SYS	66 mm ID universal fit, straight transition	3 thumbscrews
OLYMPUS	MI-150RD	MI-150 w/Ringlight ASZ Tip	fits Olympus SZ30, 40, 50	lock ring
LEICA	MI-150RG	MI-150 115V W/LEICA RING LIGHT	fits Leica S4E, S6D/E, S8APO	1 thumbscrew
UNIVERSAL	MI-150RL	MI-150 115V RINGLIGHT SYS DJI	60 mm ID universal fit, Plastic	3 thumbscrews
UNITRON	MI-150RU	MI-150 115V W/UNITRON RING LIGHT	fits Unitron ZST, ZSB (w/o aux lens)	lock ring