

## LAMP REPLACEMENT GUIDE

Model	Lamp	Wattage	Voltage	Dolan-Jenner Lamp PN		Optional Lamp	Wattage	Voltage	Lamp PN
170D	EKE	150	21	020158500000		EJV	150	21	020158520000
180	EKE	150	21	020158500000		EJV	150	21	020158520000
190	EKZ	30	10.8	020158510000					
3100	EKZ	30	10.8	020158510000					
PL750	EKE	150	21	020158500000		EJV	150	21	020158520000
PL800	EKE	150	21	020158500000		EJV	150	21	020158520000
PL900	EKE	150	21	020158500000		EJV	150	21	020158520000
DC950	EKE	150	21	020158500000		EJV	150	21	020158520000
A240L	EKE	150	21	020158500000					
A241L	EKE	150	21	020158500000					
A3000	DDL	150	20	020020000000					
A3200	DDL	150	20	020020000000					
Mi-30	EKZ	30	10.8	020158510000					
Mi-150	EKE	150	21	020158500000		EJV	150	21	020158520000
Mi-152	EKE	150	21	020158500000		EJV	150	21	020158520000
Mi-157	EKE-X	150	21	020158490000		EKE	150	21	020158500000
MH100	MHR-100	100	95	68600903162					
MVIND-150	EKE	150	21	020158500000					
I-260	MVS-2	-	n/a	020020210000					
FR-7200	LF72	15	n/a	002830030010					
FR7200B	LF72B	15	n/a	002830030020					